

Conceal

## AMS CONCEAL™

Finally, one size does fit all.

The low-profile choice is the highly flexible solution.

Conceal is the revolutionary new IPP reservoir from AMS. The low profile reservoir shape is designed to optimize fit in the Space of Retzius, offering a sleek and flexible solution for you and your patients.


- Compatible with all AMS 700 IPP cylinder (CX/LGX/CXR) and pump (700/Tactile/MS) configurations
- One size can be filled with 65-100 ml for increased flexibility
- Parylene coating technology
- Available with or without InhibiZone® antibiotic treatment

The AMS 700® with InhibiZone is the only inflatable penile prosthesis with clinical evidence showing a significant reduction in the rate of revisions due to infection in both original and revision implants.\*

\*FDA PDP D970012/S065



### RESERVOIR COMPARISON



### InhibiZone®


InhibiZone is a proprietary combination of Rifampin and Minocycline and a manufacturing process that impregnates the antibiotics into the prosthetic surfaces that contact tissue.

#### Part Numbers

720185-01 InhibiZone Conceal
720182-01 Non-InhibiZone Conceal

### InhibiZone Provides Largest Zone of Inhibition

Significantly larger zone of inhibition against MRSA than Titan,® both in vitro and in vivo, up to 14 days after implantation.¹



Zones of inhibition against Staphylococcus aureus produced by cylinder device segments explanted from rabbits. M/R = minocycline and rifampin; Van = vancomycin.

# Another innovation from AMS, a leader in men's health solutions.

1 Mansouri MD, Boone TB, Darouiche RO. Comparative assessment of antimicrobial activities of antibiotic-treated penile prostheses. Eur Urol. 56 (2009) 1039-1046.

As with any medical procedure, there are risks involved and not all patients are candidates for a penile implant. Discuss the risks and benefits of this procedure in more detail with your doctor. For a complete list of indications, contraindications and precautions, contact your urologist or refer to the Instructions for Use at www.AMSelabeling.com

